

STUKO H

Priming agent based on polymeric micro spheres, PU and acrylates for defect coverage and upgradation on leather.

PROPERTIES:

Appearance	: White Thick Paste
Base	: Microspheres
Solid Content	: 25 ± 1%
pH (10% Sol.)	: 9.0 ± 0.5
Shelf Life	: Six months

STUKO H is primer in paste form for up gradation of leather.

STUKO H is generally applied on leather with the help of blunt knife spatula and dried completely.

STUKO H treated leathers are hair cell plated and snuffed before proceeding further for finishing, thus effectively filling/hiding the surface defects of leather.

SUGGESTED APPLICATIONS:

Apply **STUKO H** with the help of blunt knife or spatula on defected parts of grain to fill up, dry properly, leave for 10-12 hours. Plate at 90°C and snuff. Proceed with the normal finishing.

FILLERS/
WAXES/
MODIFIERS

Note : Suggested formulations are only for guidance and necessary modifications must be made to achieve a particular result.